



crf1 Antibody

Product Code	CSB-PA818279ZA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8J0P4
Immunogen	Recombinant Neosartorya fumigata crf1 protein
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable glycosidase crf1 (EC 3.2.-.-) (Crh-like protein 1) (allergen Asp f 9), crf1
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	crf1
Product Modify	Unmodified