



Mcpt1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA357723ZA01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P09650 |
| Immunogen | Recombinant Rattus norvegicus Mcpt1 protein |
| Raised In | Rabbit |
| Species Reactivity | Rattus norvegicus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mast cell protease 1 (rMCP-1) (EC 3.4.21.39) (Chymase) (Chymotrypsin-like protease) (CLIP protein) (Mast cell protease I) (rMCP-I), Mcpt1 |
| Immunogen Species | Rattus norvegicus (Rat) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | Mcpt1 |
| Product Modify | Unmodified |