




Lysenin-related protein 2 Antibody

Product Code	CSB-PA517722ZA01EJV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18425
Immunogen	Recombinant Eisenia fetida Lysenin-related protein 2 protein (1-300aa)
Species Reactivity	Eisenia fetida
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Eisenia fetida
Target Names	N/A

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Lysenin-related protein 2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41.6 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.