



# G Antibody

Product Code	CSB-PA314306ZA01VBH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2X0
Immunogen	Recombinant Vesicular stomatitis Indiana virus G protein
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycoprotein G (Fragment), G
Immunogen Species	Vesicular stomatitis Indiana virus (strain Mudd-Summers) (VSIV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	G
Product Modify	Unmodified