



pep2 Antibody

Product Code	CSB-PA523645ZA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O42630
Immunogen	Recombinant Neosartorya fumigata pep2 protein
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protease A (EC 3.4.23.25) (Aspartic endopeptidase pep2) (Aspartic protease pep2), pep2
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	pep2
Product Modify	Unmodified