



map-1 Antibody

Product Code	CSB-PA765906ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WK20
Immunogen	Recombinant Mycobacterium tuberculosis map-1 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine aminopeptidase 1 (MAP 1) (MetAP 1) (EC 3.4.11.18) (Peptidase M), map-1, mapA
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	map-1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.