




malE Antibody

Product Code	CSB-PA360183ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AEX9
Immunogen	Recombinant Escherichia coli Maltose-binding periplasmic protein (malE) (27-396aa)
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli
Target Names	malE

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: malE antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 48.2 kDa

Observed band size: 49 kDa