



BLVRB Antibody, HRP conjugated

Product Code	CSB-PA002722HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30043
Immunogen	Recombinant Human Flavin reductase (NADPH) protein (2-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flavin reductase (NADPH) (FR) (EC 1.5.1.30) (Biliverdin reductase B) (BVR-B) (EC 1.3.1.24) (Biliverdin-IX beta-reductase) (Green heme-binding protein) (GHBP) (NADPH-dependent diaphorase) (NADPH-flavin reductase) (FLR), BLVRB, FLR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BLVRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.