



COL1A1 Antibody, Biotin conjugated

Product Code	CSB-PA005727OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02452
Immunogen	Peptide sequence from Human Collagen alpha-1(I) chain protein (1306-1323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(I) chain (Alpha-1 type I collagen), COL1A1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	COL1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.