



DUSP11 Antibody, Biotin conjugated

Product Code	CSB-PA007240LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75319
Immunogen	Recombinant Human RNA/RNP complex-1-interacting phosphatase protein (105-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA/RNP complex-1-interacting phosphatase (EC 3.1.3.-) (Dual specificity protein phosphatase 11) (Phosphatase that interacts with RNA/RNP complex 1), DUSP11, PIR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.