



EIF2AK2 Antibody, HRP conjugated

Product Code	CSB-PA007511LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19525
Immunogen	Recombinant Human Interferon-induced, double-stranded RNA-activated protein kinase protein (90-358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced, double-stranded RNA-activated protein kinase (EC 2.7.11.1) (Eukaryotic translation initiation factor 2-alpha kinase 2) (eIF-2A protein kinase 2) (Interferon-inducible RNA-dependent protein kinase) (P1/eIF-2A protein kinase) (Protein kinase RNA-activated) (PKR) (Protein kinase R) (Tyrosine-protein kinase EIF2AK2) (EC 2.7.10.2) (p68 kinase), EIF2AK2, PKR PRKR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EIF2AK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.