



EXOC7 Antibody, Biotin conjugated

Product Code	CSB-PA007886LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPT5
Immunogen	Recombinant Human Exocyst complex component 7 protein (521-735AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 7 (Exocyst complex component Exo70), EXOC7, EXO70 KIAA1067
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXOC7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.