



EXOC7 Antibody, Biotin conjugated

Product Code	CSB-PA007886LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPT5
Immunogen	Recombinant Human Exocyst complex component 7 protein (521-735AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 7 (Exocyst complex component Exo70), EXOC7, EXO70 KIAA1067
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXOC7