



FGF4 Antibody, HRP conjugated

Product Code	CSB-PA008631OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08620
Immunogen	Peptide sequence from Human Fibroblast growth factor 4 protein (182-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 4 (FGF-4) (Heparin secretory-transforming protein 1) (HST) (HST-1) (HSTF-1) (Heparin-binding growth factor 4) (HBGF-4) (Transforming protein KS3), FGF4, HST HSTF1 KS3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGF4