



FLT3 Antibody, HRP conjugated

Product Code	CSB-PA008733LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36888
Immunogen	Recombinant Human Receptor-type tyrosine-protein kinase FLT3 protein (727-993AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein kinase FLT3 (EC 2.7.10.1) (FL cytokine receptor) (Fetal liver kinase-2) (FLK-2) (Fms-like tyrosine kinase 3) (FLT-3) (Stem cell tyrosine kinase 1) (STK-1) (CD antigen CD135), FLT3, CD135 FLK2 STK1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FLT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.