



FZD7 Antibody, HRP conjugated

Product Code	CSB-PA009110LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75084
Immunogen	Recombinant Human Frizzled-7 protein (161-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frizzled-7 (Fz-7) (hFz7) (FzE3), FZD7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FZD7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.