



GH1 Antibody, Biotin conjugated

Product Code	CSB-PA009407LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01241
Immunogen	Recombinant Human Somatotropin protein (27-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Somatotropin (Growth hormone) (GH) (GH-N) (Growth hormone 1) (Pituitary growth hormone), GH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.