



GLI2 Antibody, HRP conjugated

Product Code	CSB-PA009500OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10070
Immunogen	Peptide sequence from Human Zinc finger protein GLI2 protein (2-20AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein GLI2 (GLI family zinc finger protein 2) (Tax helper protein), GLI2, THP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GLI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.