



HSPA1L Antibody, HRP conjugated

Product Code	CSB-PA010823HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34931
Immunogen	Recombinant Human Heat shock 70 kDa protein 1-like protein (420-641AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 1-like (Heat shock 70 kDa protein 1L) (Heat shock 70 kDa protein 1-Hom) (HSP70-Hom), HSPA1L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.