



ITGA8 Antibody, Biotin conjugated

Product Code	CSB-PA011871ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53708
Immunogen	Peptide sequence from Human Integrin alpha-8 protein (114-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-8 [Cleaved into: Integrin alpha-8 heavy chain; Integrin alpha-8 light chain], ITGA8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.