



KRT20 Antibody, Biotin conjugated

Product Code	CSB-PA012536LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35900
Immunogen	Recombinant Human Keratin, type I cytoskeletal 20 protein (207-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 20 (Cytokeratin-20) (CK-20) (Keratin-20) (K20) (Protein IT), KRT20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.