




LRPAP1 Antibody

Product Code	CSB-PA013104HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30533
Immunogen	Recombinant Human Alpha-2-macroglobulin receptor-associated protein (35-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2-macroglobulin receptor-associated protein (Alpha-2-MRAP) (Low density lipoprotein receptor-related protein-associated protein 1) (RAP), LRPAP1, A2MRAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LRPAP1

Image



Western Blot

Positive WB detected in: THP-1 whole cell lysate, A375 whole cell lysate

All lanes: LRPAP1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.