



LRPAP1 Antibody, Biotin conjugated

Product Code	CSB-PA013104HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30533
Immunogen	Recombinant Human Alpha-2-macroglobulin receptor-associated protein (35-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2-macroglobulin receptor-associated protein (Alpha-2-MRAP) (Low density lipoprotein receptor-related protein-associated protein 1) (RAP), LRPAP1, A2MRAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LRPAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.