



MFNG Antibody, HRP conjugated

Product Code	CSB-PA013757LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00587
Immunogen	Recombinant Human Beta-1,3-N-acetylglucosaminyltransferase manic fringe protein (28-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1,3-N-acetylglucosaminyltransferase manic fringe (EC 2.4.1.222) (O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase), MFNG
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	MFNG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.