



MMP2 Antibody, Biotin conjugated

Product Code	CSB-PA014666LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08253
Immunogen	Recombinant Human 72 kDa type IV collagenase protein (82-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) (TBE-1) [Cleaved into: PEX], MMP2, CLG4A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.