



MYH6 Antibody, Biotin conjugated

Product Code	CSB-PA015299HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13533
Immunogen	Recombinant Human Myosin-6 protein (894-1132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-6 (Myosin heavy chain 6) (Myosin heavy chain, cardiac muscle alpha isoform) (MyHC-alpha), MYH6, MYHCA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.