



PLA2G4B Antibody, FITC conjugated

Product Code	CSB-PA018098LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C869
Immunogen	Recombinant Human Cytosolic phospholipase A2 beta protein (514-781AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic phospholipase A2 beta (cPLA2-beta) (EC 3.1.1.4) (Phospholipase A2 group IVB), PLA2G4B
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PLA2G4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.