



RABGGTB Antibody, FITC conjugated

Product Code	CSB-PA019235HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53611
Immunogen	Recombinant Human Geranylgeranyl transferase type-2 subunit beta protein (2-331AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Geranylgeranyl transferase type-2 subunit beta (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit beta) (GGTase-II-beta) (Rab geranylgeranyltransferase subunit beta) (Rab GG transferase beta) (Rab GGTase beta) (Rab geranylgeranyltransferase subunit beta) (Type II protein geranylgeranyltransferase subunit beta), RABGGTB, GGTB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RABGGTB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.