



SERPINF1 Antibody, FITC conjugated

Product Code	CSB-PA021084LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36955
Immunogen	Recombinant Human Pigment epithelium-derived factor protein (128-418AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pigment epithelium-derived factor (PEDF) (Cell proliferation-inducing gene 35 protein) (EPC-1) (Serpin F1), SERPINF1, PEDF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.