



SESN2 Antibody, HRP conjugated

Product Code	CSB-PA021096NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58004
Immunogen	Peptide sequence from Human Sestrin-2 protein (42-61AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sestrin-2 (EC 1.11.1.15) (Hypoxia-induced gene), SESN2, Hi95 SEST2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SESN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.