



SLC3A2 Antibody, Biotin conjugated

Product Code	CSB-PA021640LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08195
Immunogen	Recombinant Human 4F2 cell-surface antigen heavy chain protein (404-630AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	4F2 cell-surface antigen heavy chain (4F2hc) (4F2 heavy chain antigen) (Lymphocyte activation antigen 4F2 large subunit) (Solute carrier family 3 member 2) (CD antigen CD98), SLC3A2, MDU1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLC3A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.