



SRPK2 Antibody, Biotin conjugated

Product Code	CSB-PA022683LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78362
Immunogen	Recombinant Human SRSF protein kinase 2 protein (471-688AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SRSF protein kinase 2 (EC 2.7.11.1) (SFRS protein kinase 2) (Serine/arginine-rich protein-specific kinase 2) (SR-protein-specific kinase 2) [Cleaved into: SRSF protein kinase 2 N-terminal; SRSF protein kinase 2 C-terminal], SRPK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SRPK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.