



TFPI2 Antibody, HRP conjugated

Product Code	CSB-PA023438LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48307
Immunogen	Recombinant Human Tissue factor pathway inhibitor 2 protein (23-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tissue factor pathway inhibitor 2 (TFPI-2) (Placental protein 5) (PP5), TFPI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TFPI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.