



TSKU Antibody, HRP conjugated

Product Code	CSB-PA025140LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUA8
Immunogen	Recombinant Human Tsukushin protein (77-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tsukushin (Tsukushi) (E2-induced gene 4 protein) (Leucine-rich repeat-containing protein 54), TSKU, E2IG4 LRRC54 TSK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TSKU
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.