



VDAC1 Antibody, Biotin conjugated

Product Code	CSB-PA025821PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21796
Immunogen	Peptide sequence from Human Voltage-dependent anion-selective channel protein 1 protein (32-59AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent anion-selective channel protein 1 (VDAC-1) (hVDAC1) (Outer mitochondrial membrane protein porin 1) (Plasmalemmal porin) (Porin 31HL) (Porin 31HM), VDAC1, VDAC
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	VDAC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.