



MDM2 Antibody, Biotin conjugated

Product Code	CSB-PA11628D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00987
Immunogen	Peptide sequence from Human E3 ubiquitin-protein ligase Mdm2 protein (374-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase Mdm2 (EC 2.3.2.27) (Double minute 2 protein) (Hdm2) (Oncoprotein Mdm2) (RING-type E3 ubiquitin transferase Mdm2) (p53-binding protein Mdm2), MDM2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MDM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.