



GPC3 Antibody, HRP conjugated

Product Code	CSB-PA12749B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51654
Immunogen	Recombinant Human Glypican-3 protein (346-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glypican-3 (GTR2-2) (Intestinal protein OCI-5) (MXR7) [Cleaved into: Secreted glypican-3], GPC3, OCI5
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	GPC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.