



AGR2 Antibody, Biotin conjugated

Product Code	CSB-PA13057D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95994
Immunogen	Recombinant Human Anterior gradient protein 2 homolog protein (21-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anterior gradient protein 2 homolog (AG-2) (hAG-2) (HPC8) (Secreted cement gland protein XAG-2 homolog), AGR2, AG2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	AGR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.