




benC Antibody

Product Code	CSB-PA27287A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1GSN3
Immunogen	Recombinant Streptomyces sp BenC protein (1-86AA)
Raised In	Rabbit
Species Reactivity	Streptomyces sp
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BenC protein, benC
Immunogen Species	Streptomyces sp. A2991200
Research Area	Others
Target Names	benC

Image



Western Blot
Positive WB detected in Recombinant protein
All lanes: benC antibody at 1:2000
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 33 kDa
Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.