




nef Antibody

Product Code	CSB-PA27307A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1L4CSJ1
Immunogen	Recombinant Human immunodeficiency virus 1 Protein Nef protein (1-209AA)
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus 1
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Nef, 3'ORF, Negative factor, F-protein, C-terminal core protein, nef
Immunogen Species	Human immunodeficiency virus 1
Research Area	Others
Target Names	nef

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: nef antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 28 kDa

Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.