



HSPA1B Antibody, FITC conjugated

Product Code	CSB-PA28047C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DMV9
Immunogen	Recombinant Human Heat shock 70 kDa protein 1B protein (418-641AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 1B (Heat shock 70 kDa protein 2) (HSP70-2) (HSP70.2), HSPA1B, HSP72
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.