




GDV1 Antibody

Product Code	CSB-PA28187A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1J1H4F2
Immunogen	Recombinant Plasmodium relictum Gametocyte development protein 1, putative protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Plasmodium relictum
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gametocyte development protein 1, putative, GDV1
Immunogen Species	Plasmodium relictum
Research Area	Others
Target Names	GDV1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: GDV1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 17 kDa

Observed band size: 17 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.