



GDV1 Antibody, HRP conjugated

Product Code	CSB-PA28187B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1J1H4F2
Immunogen	Recombinant Plasmodium relictum Gametocyte development protein 1, putative protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Plasmodium relictum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gametocyte development protein 1, putative, GDV1
Immunogen Species	Plasmodium relictum
Research Area	Others
Target Names	GDV1