



GDV1 Antibody, Biotin conjugated

Product Code	CSB-PA28187D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1J1H4F2
Immunogen	Recombinant Plasmodium relictum Gametocyte development protein 1, putative protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Plasmodium relictum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gametocyte development protein 1, putative, GDV1
Immunogen Species	Plasmodium relictum
Research Area	Others
Target Names	GDV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.