



sdeA Antibody, Biotin conjugated

Product Code	CSB-PA28219D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ZTK4
Immunogen	Recombinant Legionella pneumophila subsp Ubiquitinating/deubiquitinating enzyme SdeA protein (19-268AA)
Raised In	Rabbit
Species Reactivity	Legionella pneumophila subsp
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitinating/deubiquitinating enzyme SdeA (Effector protein SdeA) [Includes: Deubiquitinase (DUB) (EC 3.4.22.-) (Deneddylase) (Deubiquitinating enzyme); Ubiquitin transferase (EC 2.3.2.-); Mono-ADP-ribosyltransferase (mART) (EC 2.4.2.31)], sdeA
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Research Area	Others
Target Names	sdeA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.