



# TRABD2B Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA406235LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A6NFA1
<b>Immunogen</b>	Recombinant Human Metalloprotease TIKI2 protein (195-405AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Metalloprotease TIKI2 (EC 3.4.-.-) (Heart, kidney and adipose-enriched transmembrane protein homolog) (TRAB domain-containing protein 2B), TRABD2B, HKAT TIKI2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	TRABD2B