



TRABD2B Antibody, Biotin conjugated

Product Code	CSB-PA406235LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NFA1
Immunogen	Recombinant Human Metalloprotease TIKI2 protein (195-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloprotease TIKI2 (EC 3.4.-.-) (Heart, kidney and adipose-enriched transmembrane protein homolog) (TRAB domain-containing protein 2B), TRABD2B, HKAT TIKI2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRABD2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.