




USP17L3 Antibody

Product Code	CSB-PA406793LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NCW0
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 protein (19-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 3) (Ubiquitin thioesterase 17-like protein 3) (Ubiquitin-specific-processing protease 17-like protein 3), USP17L3, USP17B USP17F
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP17L3

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, K562 whole cell lysate, HepG2 whole cell lysate

All lanes: USP17L3 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 60 kDa

Observed band size: 60 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.