



USP17L3 Antibody, Biotin conjugated

Product Code	CSB-PA406793LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NCW0
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 protein (19-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 3) (Ubiquitin thioesterase 17-like protein 3) (Ubiquitin-specific-processing protease 17-like protein 3), USP17L3, USP17B USP17F
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP17L3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.