




fkpA Antibody

Product Code	CSB-PA517042HA01AXW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O08437
Immunogen	Recombinant Aeromonas hydrophila FKBP-type peptidyl-prolyl cis-trans isomerase FkpA protein (21-268AA)
Raised In	Rabbit
Species Reactivity	Aeromonas hydrophila
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	FKBP-type peptidyl-prolyl cis-trans isomerase FkpA (PPIase) (EC 5.2.1.8) (Rotamase), fkpA
Immunogen Species	Aeromonas hydrophila
Research Area	Others
Target Names	fkpA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: fkpA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.